

SUPPLEMENTARY TABLE

Supplementary Table 1. Primer sequences for real-time PCR.

Gene	Forward primer	Reverse primer	Size (bp)
CPSF6	TGGTGTGGATCTGAAGCATC	CCCAGACATTTGTCCTGATTGT	160
GSK3b	GCACTGTGTAGCCGTCTG	GAGGAGGAATAAGGATGGTAGC	199
MDM2	GAATCATCGGACTCAGGTACATC	TCTGTCTCACTAATTGCTCTCCT	167
IRS1	ACAAACGCTTCTTCGTACTGC	AGTCAGCCCGCTTGTTGATG	156
JUN	ATGGTCAGGTTATACTCCTCCTC	CACATGCCACTTGATACAATCC	164
ATF2	GTCATGGTAGCGGATTGGTTA	CTTTGGGTCTGTGGAGTTGTG	124
GAPDH	TGACTTCAACAGCGACACCCA	CACCCTGTTGCTGTAGCCAAA	121